

Suggested Treatment Areas



Philos

- FINE Lines
- DEEP Lines
- SUB-Q Lines

Acmedi KOREA

(07997) #611, Hyundai 41 Tower, 293, Mokdong dong-ro, Yangcheon-gu, Seoul, Korea
T. +82.70.5158.3466 F. +82.2.6971.9764 E. sales@acmedikorea.com
www.acmedikorea.com

Philos

Superior Quality
Cross-linked Hyaluronic Acid
for Supreme Satisfaction



Acmedi KOREA

Strong Points of *Philos* HA filler



Philos

Philos Features & Specifications

- High Purity Hyaluronic Acid (Endotoxin < 0.1EU/ml)
- Highly Cross-linked HA with HAMC (Hyaluronic Acid Moisture Correction) Technology
- Highly Consistent Mono-phasic Product - Regular & Dense Molecule Structure
- Fairly Constant Rate of Viscoelasticity Comparing to Other Products. (within 3%)



	<i>Philos</i> FINE with Lidocaine	<i>Philos</i> DEEP with Lidocaine	<i>Philos</i> SUB-Q with Lidocaine
Composition	HA 24mg/ml Lidocaine 0.3%	HA 24mg/ml Lidocaine 0.3%	HA 24mg/ml Lidocaine 0.3%
Injection Level	Superficial Dermis	Mid/ Deep Dermis	Deep Dermis Subcutaneous Layer
Recommended Indications	Superficial Lines ; Crow's feet, Forehead lines, Neck Wrinkles	Mid/Deep Folds ; Nasolabial Folds, Lips, Cheeks, Chin	Deep Folds/ Augmentation ; Face Contouring, Nose, Cheeks
Syringe Volume	1.1ml	1.1ml	1.1ml
Needle Size	30G	27G	26G
Duration	12-18 months		
Expiry	24 months from manufacture date		

Acmedi KOREA